

# Using the DAQ Control Line

Author: Jeff Landgraf

## **Components of the DAQ run control**

To run DAQ with multiple sectors you need to have a front end program - either online's DAQ Server, or the DAQ control line program. In this document, the DAQ control line is described. The front end program must connect to the DAQ Run Control Handler, which handles the communication between the front end and the components of DAQ. Currently both the front end program and the handler must run on the computer *daqman*.

## **Starting the DAQ Run Control Handler**

The Run Control Handler requires a file containing the DAQ hardware configuration. This file must exist before the Run Control Handler is started. To start the handler log into *daqman* as *daq*. Then, type:

```
cd /DAQ/bin/RC
handler -conf configfilename &
```

If you do not have a configuration file, see the section below for creating configuration files.

There is currently no message logging facility provided by DAQ. Messages are currently dumped to the file "/DAQ/log/msglog". To get a running list of the messages as they appear use the unix command:

```
tail -f /DAQ/LOG/msglog
```

## **Starting the DAQ command line:**

Once the handler has been started you are ready to start the DAQ command line. Currently only one DAQ command line session should operate at one time. To start the session type:

```
cd /DAQ/bin/RC
daqrc
```

You should get the prompt "DAQ RC→. " Press *return* to get a list of all of the valid commands.

## **Taking DATA**

To take data perform the following steps:

- 1) Type "query" at the command line. The command should return "Current State: Ready". If not proceed to the section "Fixing problems"
- 2) Type "startrun" at the command line. You will be asked to supply 3 values. For each value a default will be indicated. If you wish to accept the default value, simply press *return*. If at any time you decide not to start the run type "cancel":
  - a) Enter the run number.
  - b) Enter what kind of run you want:
    - i) Pedestal Run –

- (1) Runs pedestal calculating code for all “normal” triggers
- (2) Uses “seesaw” gain setting
- (3) Pedestal subtraction is off
- (4) Bypasses cluster finding
- (5) Writes Raw data format by default
- (6) Writes Pedestal data in the last event
- ii) Gain Run -
  - (1) //Runs Cluster Finder for “normal” events
  - (2) Uses “logarithmic” gain setting
  - (3) Pedestal Subtraction is on
  - (4) Writes GAIN data by default
- iii) Physics Run –
  - (1) //Runs Cluster Finder for “normal” events
  - (2) Uses “Gain Corrected Logarithmic” gain table // currently equivalent to “logarithmic”
  - (3) Pedestal Subtraction is set on
  - (4) Writes zero suppressed data format by default
- iv) Test Run –
  - (1) This is only for use by “DAQ experts” who know how to use it.
  - (2) Run type, pedestal mode & gain mode are taken from the “Set” values for these parameters
- c) Enter a number describing the destination for the data.
- 3) DAQ is now ready to start accepting triggers. The data will be placed in a file
  - a) If the EVD\_DEST\_MAIN variable is set to “ETHERNET” the file will be: “/data/scratch/evbData/taper1.R0#####” where # is the run number
  - b) If the EVB\_DEST\_MAIN variable is set to “SCSI” the file will be stored in: /DATA3/closed
- 4) To stop the run type “stoprun”
  - a) Currently, the program will display, “Error stopping run: There are still –1 tokens in the token manager. Do you still want to stop the run?” This is normal. Answer “y”.

## ***Setting Other Parameters***

Additional run parameters are available using the “set” command . Type “set” at the command line. This command will return a list of variables, with their current values indicated. To change a value type “set *variable\_name* value”. Only change these variables if you already understand what they mean.

## ***Fixing Problems***

## ***DAQ Run Control Handler Arguments***

Most of the handler’s command line arguments are used when running multiple instances of the handler for testing purposes. They should not be necessary for running the actual experiment. Here is the full list of command line parameters:

```
-nologtofile - tells handler not to log to a file
-log filename - tells handler filename for logging
-debug # - -1 no debug, each # * 1200 msgs
-debugtofile filnam - default daq.debug
-quiet - only log CMD_LOG_DATA messages
-lport # - give the listening port number
-padd # - add # to both listen and send ports
-psuf chars - suffix for pipes and log file
```

```

-cdir          - director for config file
-conf filnam   - filename for config file
-help         - this message

```

## ***DAQ Run Control Commands***

**startrun** – Start the run

**stoprun** – Stop the run

**pauserun** – Pause the run

**resumerun** – Resume the run

**query** – Get the current state of DAQ and write to the log file

**trigger trigger\_cmd #bunches <# in a bunch>** - Send software triggers to DAQ

**set <timeout> <variable> <value>** - Set a variable or timeout

**reconfig filename** – Change the configuration filename

**sendconfig** – Sends the current configuration to DAQ components, but does not start the run

**querytokens** – Ask all of the DAQ components how many tokens they are currently processing. The value returned is the token managers value. The values for other components are displayed by the query command

**cleardaq** – Causes all of the DAQ components to flush their buffers

**clearerror <state>**– Changes the handlers understanding of all DAQ components state to the specified state. The default state is “Ready”. You need to use the numeric state value for this command.

**reboot** – Reboot all of the connected DAQ components

**stophandler** – Exit run control and stop the handler

**quit** – Exit run control, but do not stop the handler

The following two commands allow you to control which components described in the configuration file are active in the run.

**add daqpath <rb\_mask>** – Add a component to the run

**subtract** – Remove a component from the run

## ***Example DAQ RC Hardware Configuration File***

A sample configuration file is shown below:

```

begin: experiment STAR.1
# TPC Subsystem
begin: system TPC 1

```

```
begin: detector TPC
#
begin: crate tsb12
boards: in 1 2 3 4 5 6 7 8 9 10 11 12
ip:    tsb12
end: // crate tsb12
#
begin: crate evb01
ip: evb01
evb: in
end: // crate evb01
#
begin: crate gl3
ip: in gl3
gb:
end: // crate gl3
#
end: // Detector TPC
end: // System TPC
end: // Experiment STAR.1
```

This configuration file states that there are 3 DAQ entities in the system connected to an ethernet connection: the sector broker SB12, the event builder EVB01, and the global broker GL3. It also states that the sector broker SB12 controls 12 receiver boards.